

Алматы (7273)495-231	Иваново (4932)77-34-06	Магнитогорск (3519)55-03-13	Пермь (342)205-81-47	Тверь (4822)63-31-35
Ангарск (3955)60-70-56	Ижевск (3412)26-03-58	Москва (495)268-04-70	Ростов-на-Дону (863)308-18-15	Тольятти (8482)63-91-07
Архангельск (8182)63-90-72	Иркутск (395)279-98-46	Мурманск (8152)59-64-93	Рязань (4912)46-61-64	Томск (3822)98-41-53
Астрахань (8512)99-46-04	Казань (843)206-01-48	Набережные Челны (8552)20-53-41	Самара (846)206-03-16	Тула (4872)33-79-87
Барнаул (3852)73-04-60	Калининград (4012)72-03-81	Нижний Новгород (831)429-08-12	Саранск (8342)22-96-24	Тюмень (3452)66-21-18
Белгород (4722)40-23-64	Калуга (4842)92-23-67	Новокузнецк (3843)20-46-81	Санкт-Петербург (812)309-46-40	Ульяновск (8422)24-23-59
Благовещенск (4162)22-76-07	Кемерово (3842)65-04-62	Ноябрьск (3496)41-32-12	Саратов (845)249-38-78	Улан-Удэ (3012)59-97-51
Брянск (4832)59-03-52	Киров (8332)68-02-04	Новосибирск (383)227-86-73	Севастополь (8692)22-31-93	Уфа (347)229-48-12
Владивосток (423)249-28-31	Коломна (4966)23-41-49	Омск (3812)21-46-40	Симферополь (3652)67-13-56	Хабаровск (4212)92-98-04
Владикавказ (8672)28-90-48	Кострома (4942)77-07-48	Орел (4862)44-53-42	Смоленск (4812)29-41-54	Чебоксары (8352)28-53-07
Владимир (4922)49-43-18	Краснодар (861)203-40-90	Оренбург (3532)37-68-04	Сочи (862)225-72-31	Челябинск (351)202-03-61
Волгоград (844)278-03-48	Красноярск (391)204-63-61	Пенза (8412)22-31-16	Ставрополь (8652)20-65-13	Череповец (8202)49-02-64
Вологда (8172)26-41-59	Курск (4712)77-13-04	Петрозаводск (8142)55-98-37	Сургут (3462)77-98-35	Чита (3022)38-34-83
Воронеж (473)204-51-73	Курган (3522)50-90-47	Псков (8112)59-10-37	Сыктывкар (8212)25-95-17	Якутск (4112)23-90-97
Екатеринбург (343)384-55-89	Липецк (4742)52-20-81		Тамбов (4752)50-40-97	Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

www.heathrow.nt-rt.ru | | hwe@nt-rt.ru

Технические характеристики на принадлежности для безопасности компании **Heathrow**

Виды товаров: лабораторные коврики, силиконовые коврики, пакеты для автоклавов, мешки для утилизации биологических опасностей и держатели, сводные кольца для стабилизации, бутылки для промывки.

Biohazard Disposal Bags

For personal or general lab use, sizes to suit all needs

- Polypropylene
- A range of biohazard collection bags printed with biohazard warning symbol and precautions in English, Spanish and French
- Includes an indicator patch that darkens when autoclaved



Biohazard Disposal Bags

Item No.	Color	Thickness mil	Capacity	Packaging	L x W in	L x W cm
HS10320			2 L	500	8 x 12	20.3 x 30.5
HS10321	Red ■	1.6	9 L	200	14 x 19	35.6 x 48.3
HS10322			19 L	200	19 x 23	48.3 x 58.4
HS10323			40 L	200	25 x 33	63.5 x 83.8

HDPE Coated Steel Wire Bag Holder

Item No.	L x W x H in	L x W x H cm
HS10044H	3.6 x 5 x 8.4	9.1 x 12.8 x 21.3

**Holds Biohazard Disposal Bag HS10320*



Epoxy Coated
HDPE Wire Holder



Autoclavable Disposal Bags

Clear bags imprinted with blue biohazard labeling designed for safe disposal

- Constructed of 2 mil thick polypropylene, allowing for repeated autoclaving
- Front write-on panel
- 4 inner packs of 50 bags
- Strong side seams for securely holding items



Case of 200

Item No.	Color	Capacity	L x W in	L x W cm
HS1002A		14 L	12.2 x 26	31 x 66
HS1002B	Clear □	26 L	16.3 x 23.6	41.5 x 60
HS1002C		69 L	24 x 31.9	62 x 81

Workstation Lab Mat

Safely facilitate workflow processes

- Protect counter tops from hot items up to 200°C (392°F)
- Lip around the edge can hold approx. 250 mL of liquid
- Small footprint, 12 x 14.9 x 0.1 in (30.5 x 37.9 x 0.2 cm)



Gridded Surface

for tube staging, preventing movement of objects

Reversible

Side One: combination gridded and smooth surface
Side Two: smooth surface only



Item No.	Color	L x W x H in	L x W x H cm
120747	Blue ■	12 x 14.9 x 0.1	30.5 x 37.9 x 0.2
120748	Purple ■		

Silicone Lab Mat

Large working space

The lab mat can be used with items such as centrifuges, vortex mixers, hot plates, stirrers, pipettes and so much more

- Mat features designated magnetic area to keep stir bars in place
- Protect countertops from hot items up to 200°C (392°F)



Raised Edges

enable liquid retention of approximately 300 mL of fluid

Large Working Surface

accommodates benchtop equipment, weigh boats, stir bars, scales, etc.



Item No.	Side #1 Color	Side #2 Color	L x W x H in	L x W x H cm
120506	Yellow ■	Grey/Blue ■■		
120507	Purple ■	Purple/Grey ■■	14 x 23.5 x 0.1	35 x 60 x 0.2
120599	Blue ■	White/Blue □■		

Wash Bottles Color-Coded

Color-coded for easy identification

Pack of 6

The closure and stem are molded in one piece to make them leak-proof

The angled neck and tip provide a uniform stream of liquid and the tip can be cut back if required to increase the flow

- Wide bottle mouth makes filling easier
- Bottle is made from LDPE, with a PP closure
- Color-coding for user, application or location identification

Item No.	Cap/Lid Color	Capacity
120245	White <input type="checkbox"/>	500 mL
120246	Red <input type="checkbox"/>	
120247	Yellow <input type="checkbox"/>	
120248	Green <input type="checkbox"/>	
120249	Blue <input type="checkbox"/>	



Create a safer working surface with our reusable Silicone Lab Mats

- Anti-skid to keep glassware in place
- Protect counter tops from stain, spills and wear
- Available in multiple colors and sizes



Wash Bottles Safety Labeled and Self Venting

Clear and informative product safety information

The information provided on the front of each bottle is specific to the product contained and includes the following — the Chemical Abstract Service (CAS) No., National Fire Protection Association (NFPA) codes and symbols, Protective Equipment Symbol (PES), Target Organ Information and Hazard Codes

There is also the option of using other liquids within the bottles and creating the appropriate safety information on the outside of the bottle by using the self-labeling versions

The closure and stem are molded in one piece to make them leak-proof, the angled neck and tip provide a uniform stream of liquid and the tip can be cut back to increase the flow

- Wide bottle mouth makes filling easier
- Sodium Hypochlorite bottle made from natural LDPE for light protection and has a HDPE closure
- All other bottles are made from LDPE, with a PP closure
- Non-metallic vent that prevents the risk of solvent leakage
- Six named chemicals, a self-labeling** and assorted pack of bottles to choose from

Pack of 6

Item No.	Cap/Lid Color	Description	Capacity
120250	Red	Acetone	500 mL
120251	Clear	Distilled Water	
120252	Yellow	Isopropanol	
120253	White	Ethanol	
120254	Green	Methanol	
120255	Blue	Sodium Hypochlorite	
120256	White	Self-labeling**	
120257	Assorted*	Assorted*	

* Assorted Pack includes one of each: acetone, distilled water, isopropanol, ethanol, sodium hypochlorite & a self-labeling bottle

**Self-labeling bottles must only be used with chemicals that are safe and compatible with LDPE. Labeling stickers included



Lead Rings

Stabilize vessels and labware on the bench or in water baths

- Vinyl-coated lead rings will not mark or scratch benchtops and labware
- C-shape format allows some manual size variation for a customized fit



Item No.	Shape	Color	Flask Sizes mL	ID in	ID cm	OD in	OD cm	Weight g
HS8881A	C	Blue	125–500	1.9	4.2	3	7.5	215
HS8882A	O			1.9	4.8	3	7.7	210
HS8881B	C	Blue	250–1000	2	5.1	3.5	9	465
HS8882B	O			2	5.1	3.5	9	490
HS8881C	C	Blue	500–2000	2.1	5.4	3.9	9.8	595
HS8882C	O			2.2	5.7	3.9	9.8	610
HS8881D	C	Blue	1000–4000	2.6	6.6	4.4	11.3	890
HS8882D	O			2.8	7	4.5	11.5	935

*Weights may vary by +/- 3%

Алматы (7273)495-231	Иваново (4932)77-34-06	Магнитогорск (3519)55-03-13	Пермь (342)205-81-47	Тверь (4822)63-31-35
Ангарск (3955)60-70-56	Ижевск (3412)26-03-58	Москва (495)268-04-70	Ростов-на-Дону (863)308-18-15	Тольятти (8482)63-91-07
Архангельск (8182)63-90-72	Иркутск (395)279-98-46	Мурманск (8152)59-64-93	Рязань (4912)46-61-64	Томск (3822)98-41-53
Астрахань (8512)99-46-04	Казань (843)206-01-48	Набережные Челны (8552)20-53-41	Самара (846)206-03-16	Тула (4872)33-79-87
Барнаул (3852)73-04-60	Калининград (4012)72-03-81	Нижний Новгород (831)429-08-12	Саранск (8342)22-96-24	Тюмень (3452)66-21-18
Белгород (4722)40-23-64	Калуга (4842)92-23-67	Новокузнецк (3843)20-46-81	Санкт-Петербург (812)309-46-40	Ульяновск (8422)24-23-59
Благовещенск (4162)22-76-07	Кемерово (3842)65-04-62	Ноябрьск (3496)41-32-12	Саратов (845)249-38-78	Улан-Удэ (3012)59-97-51
Брянск (4832)59-03-52	Киров (8332)68-02-04	Новосибирск (383)227-86-73	Севастополь (8692)22-31-93	Уфа (347)229-48-12
Владивосток (423)249-28-31	Коломна (4966)23-41-49	Омск (3812)21-46-40	Симферополь (3652)67-13-56	Хабаровск (4212)92-98-04
Владикавказ (8672)28-90-48	Кострома (4942)77-07-48	Орел (4862)44-53-42	Смоленск (4812)29-41-54	Чебоксары (8352)28-53-07
Владимир (4922)49-43-18	Краснодар (861)203-40-90	Оренбург (3532)37-68-04	Сочи (862)225-72-31	Челябинск (351)202-03-61
Волгоград (844)278-03-48	Красноярск (391)204-63-61	Пенза (8412)22-31-16	Ставрополь (8652)20-65-13	Череповец (8202)49-02-64
Вологда (8172)26-41-59	Курск (4712)77-13-04	Петрозаводск (8142)55-98-37	Сургут (3462)77-98-35	Чита (3022)38-34-83
Воронеж (473)204-51-73	Курган (3522)50-90-47	Псков (8112)59-10-37	Сыктывкар (8212)25-95-17	Якутск (4112)23-90-97
Екатеринбург (343)384-55-89	Липецк (4742)52-20-81		Тамбов (4752)50-40-97	Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

www.heathrow.nt-rt.ru | | hwe@nt-rt.ru